Chlamydia trachomatis

IgA EIA

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For *in vitro* diagnostic use only

Chlamydia trachomatis IgA EIA

A solid-phase enzyme immunoassay for the detection of IgA antibodies to *Chlamydia trachomatis* in human serum or plasma.

CONTENTS

	Page
INTENDED USE	2
INTRODUCTION	2
PRINCIPLE OF THE TEST	2
KIT CONTENTS	2
REAGENT PREPARATION	3
MATERIALS REQUIRED BUT NOT PROVIDED	3
SAMPLE COLLECTION AND HANDLING	3
PRECAUTIONS	4
TEST PROCEDURE	4
RESULTS	5
LIMITATIONS OF THE PROCEDURE	5
PERFORMANCE CHARACTERISTICS	5
TROUBLE SHOOTING	7
REFERENCES	21
RELATED PRODUCTS	22
SYMBOLS USED	22

INTENDED USE

Ani Lab systems' Chlamydia trachomatis IgA test has been developed for the detection of species specific IgA antibodies to *Chlamydia trachomatis* in human serum or plasma.

INTRODUCTION

Chlamydia trachomatis is a leading cause of sexually transmitted infections. It is a major cause of pelvic inflammatory disease (PID), which is a common name for different kind of infections of the uterus, fallopian tubes and adjacent pelvic structures. PID has severe long-term sequelae, specifically tubal factor infertility and ectopic pregnancy. Many PID cases are asymptomatic or minimally symptomatic but may still cause permanent tubal damage (1). Therefore screening for PID causing microorganisms, specially *Chlamydia trachomatis*, helps in prevention of these long-term consequences.

Chlamydia trachomatis infection has also been associated with reactive arthritis (2).

The accurate diagnosis of *Chlamydia trachomatis* infection by cell culture, antigen detection or DNA/PCR is often complicated and expensive. The diagnosis of a passed or complicated infection is not always possible by using the above mentioned techniques.

The alternative method for clinical diagnosis is to measure antibodies to *Chlamydia trachomatis* from serum or other secretions. However, the detection of *Chlamydia trachomatis* specific antibodies in human is complicated due to the common structures of *Chlamydia trachomatis* and *Chlamydia pneumoniae* and the high prevalence of *Chlamydia pneumoniae* antibodies in population. Common structures may lead to cross reactions in serological tests.

Ani Lab systems' *Chlamydia trachomatis* species specific EIA tests are based on synthetic peptides derived from the *C. trachomatis*-specific variable domain of MOMP (major outer membrane protein). This allows the screening and diagnosing of *Chlamydia trachomatis* infections without interference of *Chlamydia pneumoniae* antibodies (3). For complete diagnosis of acute and chronic infections it is important to measure both IgG and IgA antibodies to *Chlamydia trachomatis*.

Microimmunofluorescence (MIF) has been the common method for the *C. trachomatis* serology. It has been shown that the seroprevalence rates of peptide based EIA assays are comparable to the MIF assay (4). EIA tests are, however, less laborious, less expensive and possible to use in automates.

PRINCIPLE OF THE TEST

The principle of Ani Lab systems' Chlamydia trachomatis EIA kit is based on an indirect solid-phase enzyme immunoassay with horseradish peroxidase as the marker enzyme. The assay proceeds according to the following reactions.

When present in patient serum *Chlamydia trachomatis* IgA antibodies combine with *Chlamydia trachomatis* peptides attached to the polystyrene surface of the microplate wells. Residual patient sample is removed by washing and horseradish peroxidase conjugated anti-human IgA (sheep) is added. Unbound conjugate is washed off and a colourless enzyme substrate (H_2O_2) containing the chromogen (TMB) is added. The enzyme reaction of chromogen/substrate produces a coloured end product. Enzyme-chromogen/substrate reaction is terminated with acid (H_2SO_4). The colour intensity is directly proportional to the concentration of *Chlamydia trachomatis* antibodies in a patient sample.

KIT CONTENTS

- Reagents are stored between +2°C and +8°C.
- The expiration date is printed on each component label and on the package. Do not use reagents after the expiration date.
- Avoid excessive exposure to light. This is merely a precaution. The light sensitive reagents are the conjugate and the TMB-substrate solution, the latter one is packaged in nontransparent plastic vial for protection.
- The following reagents can be interchanged between Ani Lab systems' Chlamydia pneumoniae, Chlamydia trachomatis, Mycoplasma pneumoniae and Bordetella pertussis IgG, IgA and IgM kits and

different lots. These reagents can also be ordered separately (see Product list):

- Washing solution
- Stopping solution
- TMB- Substrate solution

- 1 MICROPLATE, 12 strips x 8 wells
Coated microplate
- 2 SAMPLE DILUENT, 30 ml
Tris buffered saline with proprietary additives, a blue colouring reagent, and 0.05 % Bronidox® as preservative.
- 3a NEGATIVE CONTROL, 1 ml
Chlamydia trachomatis antibody negative human serum with 0.05 % Bronidox® as preservative.
- 3b POSITIVE CONTROL, 1 ml
Chlamydia trachomatis antibody positive human serum with 0.05 % Bronidox® as preservative.
- 4 CONJUGATE, 40 ml
Buffered salt solution with proprietary additives, a red colouring reagent, horseradish peroxidase conjugated anti-human IgA (sheep) with 0.1% N-Methylisothiazolone as preservative.
- 5 TMB-SUBSTRATE SOLUTION, ready for use, 2 x 18 ml
Citrate buffered solution of 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide with proprietary additives and 0.01 % Kathon CG as preservative.
- 6 STOPPING SOLUTION, 25 ml
0.45 M H₂SO₄
- 7 WASHING SOLUTION, 100 ml
Concentrated citrate buffered saline, with proprietary additives, and 0.05 % Bronidox® as preservative.

INCUBATION COVERS, 2 pcs
DISPOSABLE REAGENT BASINS, 6 pcs

REAGENT PREPARATION

Reagent	Preparation	Stability of opened/ diluted reagents (+2°C to +8°C)
1 Coated microplate	Ready for use	6 months *)
2 Sample diluent	Ready for use	6 months *)
3a Negative control	Ready for use	6 months *)
3b Positive control	Ready for use	6 months *)
4 Conjugate	Ready for use	6 months *)
5 TMB-Substrate solution	Ready for use	6 months *) Discard unused reagent. A deep blue colour present in the substrate solution indicates that the solution has been contaminated and must

6 Stopping solution	Ready for use	be discarded. 6 months *)
7 Washing solution concentrate (10x)	Dilute the concentrate (vial 7)	6 months *)
Washing solution	1+9 (1:10) with distilled water	1 month at +4°C or 1 week at room temperature

*) Once the microplate foil-package is opened it should be resealed tightly with the desiccant: fold the opened end a few times and seal air-tightly with tape over the whole length of the opening. The stability of the opened reagents is the maximum only if they are stored properly at +2°C to +8°C. High environmental temperature and contamination may decrease the stability.

MATERIALS REQUIRED BUT NOT PROVIDED

- Distilled or deionized water, preferably sterile.
- Graduated cylinders for reagent dilution.
- Vials to store the diluted reagents.
- Precision pipettes (one channel eg. 0.5-10 µl, 5-50 µl, 20-200 µl, 100-1000 µl ranges and multi-channel 50-300 µl)
- Paper towels or absorbent paper.
- Timer, 60 min range.
- Microplate incubator
- Microplate photometer, 450 nm
- Microplate washer
- Sodium hypochlorite solution, free available chlorine 50-500 mg/l.
- Disposable gloves.

SAMPLE COLLECTION AND HANDLING

Chlamydia trachomatis IgG and IgA EIAs can be performed only from serum and plasma (EDTA, Li-Heparin) samples. Paired samples should be collected in the similar way. The use of citrate and ACD plasma is not allowed due to dilution of the plasma with the anticoagulant which may affect the clinical interpretation.

Serum and plasma samples should be refrigerated (+4°C) after collection or, if the test cannot be performed within 1 week, frozen (-20°C or -70°C, which is preferred). **Samples should not be repeatedly frozen and thawed.**

Do not use sodium azide as preservative because it inactivates horseradish peroxidase.

Heat inactivation of serum or plasma (+56°C, 30 min.) may cause an unspecific result.

Microbially contaminated, grossly hemolyzed or hyperlipemic serum and plasma may give erroneous results. Long storage of serum (frozen over one year) may cause the formation of lipid aggregates. These aggregates may cause an unspecific result.

PRECAUTIONS

For *in vitro* diagnostic use only.

Warning - POTENTIAL BIOHAZARDOUS MATERIAL:

All human materials used in the preparation of the calibrators/controls in the kit have been tested for the presence of the antibodies to HIV (Human Immunodeficiency Virus) and HCV (Hepatitis C Virus) as well as Hepatitis B surface antigen (HBsAg) and found to be non-reactive. Because no test method can offer complete assurance that HIV, hepatitis B virus, HCV, or other infectious agents are absent, these calibrators and controls as well as samples should be handled at the Biosafety level 2 as recommended for any potentially infectious human serum or blood sample in the Centers for Disease Control/National Institutes for Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 2007 (6).

Discard all materials and samples as if capable of transmitting infection. The preferred method of disposal is autoclaving for a minimum of one hour at 121°C. Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 50-500 mg/l free available chlorine. Allow 30 minutes for decontamination to be completed. Spills should be wiped off thoroughly using either an iodophor disinfectant or sodium hypochlorite solution. Materials used to wipe off spills should be added to biohazardous waste matter for proper disposal. Reusable glassware must be disinfected, washed out and rinsed free of detergents.

Liquid waste containing acid must be neutralized with a proportional amount of base prior to the addition of sodium hypochlorite. Stopping solution (vial 6) contains 0.45 M sulphuric acid, avoid contact with skin and eyes.

Wear disposable gloves while handling samples and kit reagents. Afterwards wash hands carefully. Never pipette by mouth.

Do not mix or interchange reagents from different lots, except washing solution, stopping solution, and TMB-substrate solution. Do not interchange vial caps.

Once the assay has been started, all subsequent steps should be performed without interruption. Do not let the wells dry once the assay has been started.

Do not reuse a strip of the microplate even if some wells were not used.

Accurate and precise pipetting, as well as following the exact time and temperature requirements, is essential.

TEST PROCEDURE

Preliminary preparations

- Bring the reagents and microplates to room temperature (+20°C to +25°C) before starting the assay.
- Prewarm the incubator to +37°C.

- Note: When using an automate, please check that the real incubation times are not too long

STEP I

- Pipette 180 µl of the sample diluent (vial 2) into each well²
- Pipette 20 µl of sample diluent as reagent blank, 20 µl of samples, 20 µl of negative control (vial 3a) and 20 µl of positive control (vial 3b). Mix well.¹
- Cover the plate with plastic sheet. Incubate for 30 (± 3) minutes at +37°C (±1 °C).
- Wash 5 x 400 µl / well³

STEP II

- Pipette 200 µl of the conjugate (vial 4) into each well²
- Cover the plate with plastic sheet. Incubate for 30 (± 3) minutes at +37°C (±1 °C).
- Wash 5 x 400 µl / well³

STEP III

- Add 200 µl of the TMB-substrate solution (vial 5) into each well.²
- Incubate for 15 (±1) minutes at room temperature (+20-25 °C) in dark.

STEP IV

- Stop the enzyme-substrate reaction by adding 100 µl of stopping solution (vial 6) into each well.⁴
- Measure the absorbances immediately at 450 nm /reference 620 nm (590 -690 nm)⁵

NOTES

The use of an 8-channel pipette device is recommended for improved efficiency and precision.

1. Use of duplicates is preferable, especially for the Positive control. After pipetting of the sample diluent into the wells, pipette first the samples, and then the controls. Use a new pipette tip for each sample. Mix well: insert the pipette tip into sample diluent and mix the sample/control gently while pipetting to avoid gradient. Do not touch the walls of the wells with pipette tips when mixing the samples

2. **Avoid contamination:** When removing aliquots from the reagent vials, use aseptic technique to avoid contamination. Pour needed amount of the sample diluent, **conjugate and TMB-substrate solutions into a disposable basin for reagents** (6 pieces supplied with the kit). **Discard any unused solutions; do not pour it back to the vials.** Do not touch the walls of the wells with pipette tips when adding TMB-substrate.

3. Washing may be performed manually or with a washer. Washing solution is recommended to stay in wells for 15-30 seconds during each cycle. After the washing step tap the inverted plate a few times on the paper towel.

4. When stopping reagent is added the blue colour of positive reaction should turn yellow. Green colour indicates that the stopping reagent is not properly mixed in the well. If needed mix the wells before measuring

5. The absorbance of the reagent blank has to be measured to check that it falls within the Quality Control values

RESULTS

Quality Control Values

Before calculating the results, make sure that the absorbance values obtained for the reagent blank and controls fall within the Quality Control guidelines. If the blank and the controls do not give expected values, the results are invalid and the samples should be retested.

Quality control values

QC Sample	Expected value at 450 nm (in absorbance units)
Reagent blank	< 0.10
Negative control NC	< 0.20 *)
Positive control PC	0.70 ≤ PC ≤ 2.00 *)

*) = The reagent blank absorbance has already been subtracted from these values.

Calculation of the Results

Appreviations:

- Arb = Mean absorbance of the reagent blank
- Apc = Mean absorbance of the positive control
- As = Mean absorbance of the sample
- CO = Cut-off value in absorbance units

Use the following formula for calculating the **CO** (cut-off) and **S/CO** (signal to cut-off) values:

$$CO = 0.6 \times (Apc - Arb)$$

$$S/CO = (As - Arb) / CO$$

Example:

Sample	Mean A at 450 nm
Reagent blank	0.030
Positive control	1.370
Sample 1	1.520

$$CO = 0.6 \times (1.370 - 0.030) = 0.804$$

$$S/CO \text{ (Sample 1)} = (1.520 - 0.030) / 0.804 = 1.85$$

Interpretation of the Results

- S/CO < 1.0** Negative.
A negative result indicates that the sample tested either contains no IgA antibodies to *Chlamydia trachomatis* or that the antibody level is not yet detectable. Sample may contain IgA antibodies to *Chlamydia trachomatis*.
- 1.0 ≤ S/CO < 1.4** Equivocal.
An equivocal result indicates that presence of anti-*Chlamydia trachomatis* IgA antibodies can not be definitely established. In case of clinical suspicion retest antibodies after two weeks.
- S/CO ≥ 1.4** Positive.
A positive result may indicate either past or acute on-going infection.

LIMITATIONS OF THE PROCEDURE

Because no single method leads to the definitive diagnosis, the results of the present method should be interpreted in conjunction with the clinical condition and other laboratory methods.

It is recommended that the assay is performed by qualified and trained laboratory technician.

PERFORMANCE CHARACTERISTICS

Reproducibility

Table 1. Within run reproducibility

C. trachomatis IgA Sample no	replicates n	mean absorbance	SD	CV %
1	24	1.051	0.073	6.9
2	24	1.465	0.085	5.8
3	24	1.607	0.106	6.6

Table 2. Between run reproducibility

C. trachomatis IgA Sample no	replicates within-run	run n	mean abs	mean S/CO	SD S/CO	CV %
1	4	10	1.077	1.12	0.04	3.7
2	4	10	1.204	1.25	0.07	5.7
3	4	10	1.550	1.62	0.10	6.4
4	4	10	1.909	1.99	0.16	7.7

Summary of the evaluation studies

Diagnostic sensitivity

The diagnostic sensitivity of the Ani Lab systems' species specific *Chlamydia trachomatis* IgG and IgA EIAs was compared to another commercial *Chlamydia trachomatis* IgG and IgA assay using serum samples (n=56) from clinically verified Pelvic Inflammatory Disease (PID) patients.


Figure 1. Correlation of S/CO values between the Ani Lab systems C. trachomatis IgG EIA and another peptide-based method.

Table 3. The agreement of interpretation between the Ani Labsystems C. trachomatis IgG EIA and another peptide-based method.

IgG		Labsystems	
		Pos	Neg
Comparative	Pos	24	2
	Neg	6	24


Figure 2. Correlation of S/CO values between the Ani Labsystems C. trachomatis IgA EIA and another peptide-based method.

Table 4. The agreement of interpretation between the Ani Labsystems C. trachomatis IgA EIA and another peptide-based method.

IgA		Labsystems	
		Pos	Neg
Comparative	Pos	9	2
	Neg	11	34

The diagnostic sensitivities of Ani Labsystems' Chlamydia trachomatis IgG and IgA EIAs together with another commercial peptide-based C.trachomatis IgG and IgA EIA assays were evaluated in a study by Verkooyen and his group (5). Sera from patients with PCR-proven *Chlamydia trachomatis* infection were tested as well as a group of blood donors as a negative control group.

In a group of patients with PCR-proven *Chlamydia trachomatis* infection about 70 % showed positive IgG and less than 50 % positive IgA response (Table 5).

In this study less than 100 % of PCR-positive samples had *C. trachomatis* antibodies. A possible explanation for this is that the immune system has not had enough time to response.

Table 5. The prevalence of *Chlamydia trachomatis* IgG and IgA antibodies using peptide-based EIA assays from Ani Labsystems and from another commercial manufacturer (5)

	Patients with PCR-confirmed active <i>trachomatis</i> infection, n=324	Blood donors, n=443
	Ctra IgG positive (%)	
Ani Labsystems EIA	69	6
Comparative EIA	68	6
	Ctra IgA positive (%)	
Ani Labsystems EIA	38	5
Comparative EIA	48	5

Diagnostic specificity

The diagnostic specificity of the Ani Labsystems species specific Chlamydia trachomatis IgG and IgA EIAs was compared to another commercial C. trachomatis IgG and IgA assay using serum samples from children 1-5 years old (n=123). The respective results are in Fig. 3 and 4.


Figure 3. Correlation of S/CO values between the Ani Labsystems C. trachomatis IgG EIA and another peptide-based method.


Figure 4. Correlation of S/CO values between the Ani Lab systems C. trachomatis IgA EIA and another peptide-based method.

Conclusion

The Ani Lab systems C. trachomatis IgG and IgA kits show good diagnostic sensitivity and specificity.

TROUBLE SHOOTING

Cause/Error	Solution
BLANK IS TOO HIGH	
1. Contamination, spills from other wells	Avoid contamination
2. Washing solution concentrate was not diluted correctly	Should be diluted 1:10 (1+9)
3. Poor washing	Check your washer
4. Contamination of the TMB-substrate in the reaction basin.	Keep the residual substrate solution until the test is completed. Check if the substrate in the reagent basin turns blue. This will indicate contamination.

Cause/Error	Solution
ABSORBANCE VALUES ARE LOW	
1. Incubation temperature is too low	Incubate at 37°C ±1°C. The heating efficiency of different incubators vary widely. The incubator with efficient and homogeneous warming is preferred
2. Contamination of conjugate with spills of human sera.	All absorbance values are low. Only nanoliters of sera are enough to block the activity of the conjugate. Pay special attention to prevent contamination
3. Contamination of wells with spills of human sera.	Few individual absorbance values are (very) low. Pay special attention to prevent contamination

4. Reagents are deteriorated *due to contamination or HRP deterioration in the conjugate *due to improper storage	When removing aliquots from the reagent vials, use aseptic technique to avoid contamination or erroneous results may occur. Protect reagents from excessive light Store at +4°C
5. Reagents are not warmed up to room temperature before starting	Should be +20...+25°C when starting the assay
6. Incubation time is too short	Incubate 1h (+/-5min) or 30 min.
7. Interchange of reagents	Do not mix or interchange reagents from different lots
8. Stopping solution is not mixed properly	Mix carefully before measurement

POSITIVE CONTROL HAS TOO LOW ABSORBANCE VALUES	
Cause/Error	Solution
1. The control is not mixed with the sample diluent while pipetting	After dispensing the control rinse the pipette tip with sample diluent

Cause/Error	Solution
POOR PRECISION	
1. Pipettes are not properly calibrated	Check calibration of the pipettes
2. Improper washing due to contamination of tips of the washer	Clean regularly tips of the washer
3. The plate is allowed to stay too long after washing (drying of the plate)	Follow strictly the kit instructions
4. Heterogeneous warming of the plate	Service Incubator /Shaker or incubator
5. Sample serum/plasma is not mixed properly with sample buffer	While pipetting mix the sample with sample buffer
6. Stopping solution has not been mixed properly before measurement	Mix the plate before measuring

TOO HIGH ABSORBANCE VALUES	
Cause/Error	Solution
1. Incubation temperature is too high	Incubate at 37°C (± 1°C). The heating efficiency of incubators varies widely. The incubator with efficient and even warming is preferred
2. Incubation time is too long. This may be the case especially with automatic systems.	Check the real incubation times (should be 30 (±3) min, 30 (±3) min, 15 (±1) min)
3. Negative control is too high. Contamination, spills from other wells, interchange of the vial caps	Avoid contamination

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REFERENCES / LITERATUR /LITERATUR

1. Paavonen J.(1998). Pelvic inflammatory disease. From diagnosis to prevention. *Dermatologic Clinics* 16:747-756
2. Nikkari S, Puolakkainen M, Närvänen A, Aakre O, Toivanen P, Leirisalo-Repo M (2001). Use of Peptide Based Enzyme Immunoassay in Diagnosis of *Chlamydia trachomatis* Triggered Reactive Arthritis. *J Rheumatol* 28:2487-2493 .
3. Närvänen A, Puolakkainen M, Wu H, Kino K, Suni J. (1997). Detection of Antibodies to *Chlamydia trachomatis* with Peptide-Based Species-Specific Enzyme Immunoassay. *Inf. Dis. Obst. Gyn.* 5: 349-354.
4. Morré SA, Munk C, Persson K, Krüger-Kjaer S, van Dijk R, Meijer CJLM, van den Brule AJC. (2002). Comparison of three commercially available peptide-based immunoglobulin G (IgG) and IgA assays to microimmunofluorescence assay for detection of *Chlamydia trachomatis* antibodies . *J.Clin. Microbiol.* 40:584-587.
5. Verkooyen RP, Peeters MF, van Rijsoort-Vos JH, van der Meijden WI, Mouton JW (2002). Sensitivity and specificity of three new commercially available *Chlamydia trachomatis* tests. *Int. J. STD AIDS* 13 (Suppl.2)23-25.
6. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention and National Institutes of Health. 5th Edition 2007. US Government Printing Office. Washington 2007

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